


- |                                                                                             |                                                                                    |
|---------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|
| <b>(D)</b> Gebrauchsanweisung/Technische Beschreibung<br>ACCULAN-Reamer GA 613              | <b>(I)</b> Istruzioni per l'uso/Descrizione tecnica<br>Alesatore ACCULAN GA 613    |
| <b>(GB)</b> Instructions for use/Technical description<br>ACCULAN reamer GA 613             | <b>(P)</b> Instruções de utilização/Descrição técnica<br>Escareador ACCULAN GA 613 |
| <b>(F)</b> Mode d'emploi/Description technique<br>Moteur de fraisage/alésage ACCULAN GA 613 | <b>(NL)</b> Gebruiksaanwijzing/Technische beschrijving<br>ACCULAN ruimer GA 613    |
| <b>(E)</b> Instrucciones de manejo/Descripción técnica<br>Escariador ACCULAN GA 613         | <b>(S)</b> Bruksanvisning/Teknisk beskrivning<br>ACCULAN mörghorr GA 613           |

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Read these instructions for use carefully before putting the instrument into operation. This will enable you to avoid damage caused by improper assembly or operation and therefore not covered by the guarantee and warranty.

The ACCULAN reamer and the accessories may be used only by qualified personnel familiar with the surgical techniques.

We reserve the right to change design and technical data without prior notice.

Combinations of accessories other than those stated in the instructions for use may be used only if they do not impair the features and safety requirements of the intended applications.

Only Aesculap products may be combined!

## 1. Overview

The operational ACCULAN reamer GA 613 consists of a charged battery 9, the cover 7 and the reamer 1.

The integrated HARRIS coupling of the cutter permits the use of Aesculap tools with a HARRIS connector.

In addition, four further longitudinally drilled adapters permit universal use:

- Hudson GB 421
- AO-Protek GB 422
- DIN 58809 (medullary reamer) GB 423
- Self-centring chuck GB 424  
(Clamping range 0.5 – 8 mm)

- 1 ACCULAN reamer
- 2 Clamping sleeve for clamping the tools
- 3 Trigger for speed control
- 4 Gear changer
- 5 Release button for gear changer
- 6 Locking pin
- 7 Cover GA 624
- 8 Cover release
- 9 Battery GA 626
- 10 Funnel GA 625
- 11 Charging side
- 12 Drive surface of the HARRIS shaft

### Note

The ACCULAN reamer GA 613 must be operated only with original Aesculap accessories.

## Explanation of symbols

|  |                                                                                                                                                                                                                          |
|--|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Attention, see Instructions for use                                                                                                                                                                                      |
|  | Classification type B                                                                                                                                                                                                    |
|  | CE marking according to Directive 93/42/EEC                                                                                                                                                                              |
|  | Protect from moisture                                                                                                                                                                                                    |
|  | Speed control<br>Trigger not depressed: (0)<br>Machine not operating (0)<br>The speed increases as the trigger is depressed to an increasing extent<br>Trigger fully depressed: (I)<br>Machine runs at maximum speed (I) |
|  | 1st gear acetabulum cutting                                                                                                                                                                                              |
|  | 2nd gear medullary reaming                                                                                                                                                                                               |
|  | Marking of electric and electronic devices according to directive 2002/96/EC (WEEE), see Disposal.                                                                                                                       |

## Contents

|                                                  |    |
|--------------------------------------------------|----|
| 1. Overview.....                                 | 8  |
| 2. Purpose.....                                  | 9  |
| 3. Mode of operation.....                        | 9  |
| 4. Putting into operation .....                  | 9  |
| 5. Installation.....                             | 9  |
| 5.1 Inserting the battery (Fig. 1 – 3) .....     | 9  |
| 5.2 Removing the battery (Fig. 4).....           | 9  |
| 5.3 Operation with Cement Injector GB 171 N..... | 9  |
| 6. Operation .....                               | 10 |
| 7. Function test .....                           | 10 |
| 8. Care and handling .....                       | 10 |
| 8.1 Disinfecting and cleaning .....              | 10 |
| 8.2 Care and checking for damage.....            | 11 |
| 8.3 Sterilization method and parameters .....    | 11 |
| 9. Maintenance .....                             | 11 |
| 10. Instructions for use of battery GA 626 ..... | 12 |
| 11. Troubleshooting list.....                    | 12 |
| 12. Technical Service.....                       | 13 |
| 13. Technical data .....                         | 13 |
| 14. Accessories.....                             | 13 |
| 15. Disposal .....                               | 13 |

## 2. Purpose

The ACCULAN reamer GA 613 is used for cutting an acetabulum and for medullary reaming in bone surgery.

## 3. Mode of operation

In the ACCULAN reamer, the electric motor is supplied with power by a battery and the output speed is electronically controlled. The reamer is equipped with 2 gears. In 1st gear, the speed can be regulated up to 180 rpm with the trigger. In 2nd gear, the reamer reaches a speed of max. 380 rpm.

The reamer has a right-handed drive.

The reamer is provided with a hole for receiving guide spikes up to a diameter of 3.8 mm.

### Note

*The continuous speed regulation necessitates electronic components. The conditions prevailing during sterilization (steam, temperatures, etc.) subject these components to considerable stress. To avoid this stress, the entire electronics of the ACCULAN reamer are integrated in the battery 9, which must not be sterilized. The battery 9, too, should therefore be treated with the care usual for a surgical instrument.*

## 4. Putting into operation

To put the ACCULAN reamer GA 613 into operation, you require the following components:

- ACCULAN charger GA 627 M
- Battery GA 626
- ACCULAN reamer GA 613

Before putting the reamer into operation, ensure that the required accessories are available. You require charger GA 627 M and at least one battery GA 626 to guarantee operation of the ACCULAN reamer.

### Note

*The special Aesculap Instructions for use must be followed for putting the ACCULAN charger into operation.*

*Please also see the separate instructions for handling and care of the battery GA 626 under Section 10.*

## 5. Installation

### 5.1 Inserting the battery (Fig. 1 – 3)

#### ACHTUNG

Nickel-cadmium batteries must not be heated above 65 °C and therefore must not be sterilized.

The following procedure should be adopted for the sterile preparation of the reamer:

The reamer is placed upside down (battery shaft pointing upwards) and the funnel 10 (sterile) is mounted on the reamer. The battery 9 (nonsterile) is then inserted into the charging shaft by a second person (the battery cannot rotate) and the funnel 10 is removed and is placed in the nonsterile room. The cover 7 (sterile) is then mounted and must be snapped in securely.

The cover provides an absolute seal so that there is no danger of contamination of the operation site by the unsterile battery.

### 5.2 Removing the battery (Fig. 4)

After the end of the operation and prior to cleaning and disinfection and especially before sterilization, the battery must be removed.

To do this, the two cover releases 8 on the cover 7 must be pressed simultaneously to enable the cover to be removed (the cover cannot be opened if only one release is pressed). This ensures that accidental opening is prevented if only one release is pressed. The battery can be removed from the battery shaft by tapping on the palm of the hand (Fig. 4).

Prior to battery removal, make sure that the trigger for speed control 3 is in the front position (rest position).

#### IMPORTANT

Knocking against hard objects damages the machine and may damage the battery too.

### 5.3 Operation with Cement Injector GB 171 N

#### CAUTION

The battery GA 626 must always be changed over after cutting or reaming using the ACCULAN reamer GA 613 and just prior to using the cement injector. This is because only a fully charged battery is sufficient to ensure that a change of battery during cement injection can be avoided.

#### Note

*Removal of the battery GA 626 intraoperatively should be prepared using the battery changer GA 614. Please refer to section 5.1 in the users instructions.*

## 6. Operation

### Prevention of unintentional operation



In order to prevent unintentional operation (for example during tool change), the machine is equipped with a lock on the speed control trigger.

### Explanation of symbols:

ON: The machine is in the ready state.

OFF: The speed control trigger is blocked.  
The machine cannot be started.

- a By turning the trigger lock to the OFF position, you can block the speed control trigger.
- b By turning the trigger lock to the ON position, the speed control trigger is released. The machine is now in the ready state.

By pressing the trigger, the reamer is controlled continuously in 1st gear (cutting) up to 180 rpm. In 2nd gear (medullary reaming), the machine reaches max. 380 rpm.

### Note

*Owing to the electronics used, a quiet beep tone is heard shortly before the gentle start-up of the reamer. However, this tone becomes increasingly quiet as the reamer runs up.*

Switching between the gears is performed by pressing the release button 5 and then turning the gear changer 4 through 180° until the changer snaps in detectably.

### Note

*Changing between the gears is possible both when the machine is non-operational and when it is running. If the gear changer is not locked in one of the end positions, the motor runs when the trigger is pressed but the spindle with the tool remains stationary.*

The ACCULAN reamer is equipped with a HARRIS chuck. You can couple Aesculap tools with a HARRIS shaft.

To clamp the tools, pull back the clamping sleeve 2, now insert the tool shaft and turn it gently until the drive surfaces 12 snap in. Lock the tool and the handpiece by releasing the clamping sleeve 2.

To release the tools, push back the clamping sleeve 2 and pull out the tool shaft.

## 7. Function test

Each time before use, the ACCULAN reamer with the mounted tool must be subjected to a test run. Before the function test, ensure that neither the machine nor the battery is mechanically damaged. During the test run, ensure that the tools fit correctly on the reamer. For this function test, the battery must be inserted beforehand (see 5.1).

### Note

*Each intraoperative battery change has to be followed by a function check!*

## 8. Care and handling

### 8.1 Disinfecting and cleaning

#### CAUTION

Risk of damage to the ACCULAN reamer from fluid inside the device!

- Never lay the device in fluid.
- Never spray pressurized air on the device.
- Allow any fluid that has penetrated the device to run out immediately.

- Remove the tool and/or adapter.

#### Manual cleaning and disinfecting

- To clean the surfaces of the ACCULAN reamer, brush them off lightly under running water or clean them with an alcohol solution.
- Wipe off the motor with a clean, lint-free cloth that has been soaked in cleaning and disinfectant solution.
  - or -
  - In order to avoid formation of sticky or encrusted residues, spray the components with a disinfectant spray.
- Wipe off the component surfaces with a clean, lint-free cloth.
- Follow the disinfectant manufacturer's instructions.

## Mechanical cleaning and disinfecting

### CAUTION

Improper mechanical cleaning will damage the ACCULAN reamer!

- Clean the ACCULAN reamer in a suitable storage tray only.
- Use a suitable cleaning and disinfecting cycle for mechanical cleaning.
- Do not use alkaline cleaning agents.

- Perform mechanical cleaning immediately after the operation.
- Place storage tray GB 509 R in a suitable wire basket (e.g. JF 214 R).
- Position the ACCULAN reamer properly in storage tray GB 509 R immediately after the operation.
- To avoid abrupt changes in temperature, use a single-chamber washer disinfectant.
- Select a cleaning and disinfecting cycle according to the following criteria:
  - Neutral cleaner
  - Cleaning temperature ≤ 50 °C
  - Final rinse with simultaneous thermal disinfection in demineralized water.
- Example: Miele Vario TD process.
- Follow the instructions for use of the Aesculap ECCOS storage system TA 009 721.

## 8.2 Care and checking for damage

- Each time the device is processed and a function check is carried out, check for:
  - damage, atypical sounds when the device is operated or excessive vibration or heating-up
  - damage to the battery compartment seal
  - damaged or missing silicone seal on the cover
- If any irregularities are found, do not make further use of the device.
- Prior to an extended period of non-use, always store the device cleaned and dried according to the instructions.

## 8.3 Sterilization method and parameters

The ACCULAN reamer can be sterilized in the storage tray.

### CAUTION

If the rechargeable battery and its integrated electronic component are sterilized, they will be damaged!

- Remove the rechargeable battery from the device prior to sterilization.
- Do not heat the rechargeable battery to a temperature of greater than 65 °C.
- Do not sterilize the rechargeable battery.
- Sterilize the ACCULAN reamer immediately after cleaning it.
- Removing the rechargeable battery 2 from the device, see Removing the battery (Fig. 4).
- Store the ACCULAN reamer, battery compartment cover, sterile funnel, adapter and tool in a suitable wire basket (e.g., JF 214 R).
- Sterilize using steam sterilization, taking note of the following: The sterilization has to be done according to a validated steam sterilization procedure (e.g. in a sterilizer according to EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993 and validated according to EN 554/ISO 13683). In case of application of the fractionated vacuum procedure the sterilization has to be done using the 134 °C/2 bar program with a minimal holding time of 5 minutes.

## 9. Maintenance

To ensure reliable operation, Aesculap recommends that maintenance be performed after 300 processing cycles, at least once per year.

- For service, maintenance and repairs, contact the Aesculap/B. Braun representative in your country, see Technical Service.

### Note

*The battery compartment of the reamer 1 is closed germ-tight and water-tight with the cover 7. The silicone seal in the cover 7 and the sealing surface on the reamer handle should therefore be checked regularly. A damaged seal and/or sealing surface may lead to the introduction of germs from the battery shaft.*

*Aesculap's Technical Service should be consulted even in the event of slight damage. This enables you to avoid greater damage or the failure of your handpieces.*

## 10. Instructions for use of battery GA 626

The battery GA 626 is a powerful battery with integrated electronics and power source. Owing to its special design and the computer-controlled charger, the battery has a long life if it is properly treated.

### IMPORTANT

Avoid heating the battery to temperatures above 65 °C. Heating above this temperature causes the battery capacity to decrease and prolonged exposure to this temperature destroys the battery.

### IMPORTANT

Avoid short-circuiting the battery. If the poles of the battery (top or bottom) are short-circuited by liquids or metal parts, a very high current (roughly comparable with that of a car battery) flows. The battery may therefore heat up to a very great extent and also damage adjacent objects through the action of heat (there is also a danger of fire).

We therefore advise the following handling in your own interest and the interest of your ACCULAN system:

- Always store the battery at room temperature.

- Do not sterilize or disinfect the battery (overheating, danger of short-circuit.)  
Owing to the complete encapsulation of the battery in the machine, there is absolutely no risk of contamination of the operation site by the battery. Disinfection or sterilization is therefore not necessary.
- Always keep the battery dry.
- Never immerse the battery in liquid.
- Never place the battery on a moist surface.
- Before you place the battery in the charger, ensure that the inside of the charging shaft is dry. If this is not the case, dry it with a cloth.
- Do not store or use the charger and battery in the vicinity of workplaces where large quantities of liquid are handled (instrument cleaning, transfer of wash liquid, etc.).

### Protect from moisture

If you observe these measures, you will certainly be very pleased with your ACCULAN system.

If you have any questions about the ACCULAN system, consult your medical product adviser.

## 11. Troubleshooting list

| Fault                                         | Cause                                                                                                                                                                                                                                                | Fault detection                                                             | Rectification                                                 |
|-----------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------|---------------------------------------------------------------|
| Motor runs too slowly or has too low a torque | NiCd battery is too weak<br>NiCd-Akku is defective or flat                                                                                                                                                                                           |                                                                             | Change battery<br>Check charger, if necessary replace battery |
| Motor does not run                            | If the ACCULAN reamer is used infrequently and in addition the applications are relatively brief, an oxidation film may form on the motor/battery connectors. This oxidation film is an insulating layer through which the low voltage may not pass. | No beep tone is heard from the electronics when the trigger is pressed      | Briefly turn the drive shaft by hand                          |
|                                               | Gear or bearing defective                                                                                                                                                                                                                            | The beep tone from the electronics can be heard when the trigger is pressed | Repair by manufacturer                                        |
| Motor runs but tool does not turn             | Gear changer is not snapped in at the end position                                                                                                                                                                                                   | Motor runs audibly                                                          | Snap in gear changer at end position                          |

## 12. Technical Service

For service, maintenance or repairs contact your Aesculap/B. Braun representative.

All guarantees and warranty rights as well as licenses will be void if modifications on the device are performed.

### Service addresses

Aesculap Technischer Service  
Am Aesculap-Platz  
78532 Tuttlingen / Germany  
Phone: +49 7461 95-2700  
Fax: +49 7461 16-2887  
E-mail: ats@Aesculap.de

Other service addresses can be obtained from the address indicated above.

## 13. Technical data

|                               |                         |             |
|-------------------------------|-------------------------|-------------|
| Machine type                  | GA 613                  |             |
| DC voltage                    | 9.6 V                   |             |
| Rated battery capacity        | 0.8 Ah                  |             |
| Rated current                 | 6 A                     |             |
| Max. mechanical power         | 80 W                    |             |
|                               | 1st gear                | 2nd gear    |
| Oscillation frequency         | 0 – 180 rpm             | 0 – 380 rpm |
| Max. torque                   | 15 Nm                   | 7 Nm        |
| Dimensions (H x B x T)        | 225 mm x 190 mm x 55 mm |             |
| Weight (with battery)         | 1600 g                  |             |
| Type (according to IEC 601)   | B                       |             |
| Electromagnetic compatibility | IEC 601-1-2             |             |
| Standards applied             | IEC 601-1               |             |

## 14. Accessories

|                                                  |          |
|--------------------------------------------------|----------|
| Charger                                          | GA 627 M |
| Battery                                          | GA 626   |
| Cover                                            | GA 624   |
| Sterile funnel                                   | GA 625   |
| Longitudinally drilled adapters Hudson           | GB 421   |
| AO-Protek                                        | GB 422   |
| DIN 58809 (medullary reamer)                     | GB 423   |
| Self-centring chuck<br>Clamping range 0.5 – 8 mm | GB 424   |
| Storage tray (for GA 613)                        | GB 509 R |
| Storage tray (for GA 624)                        | GB 504 R |

## 15. Disposal

### Note

*The user institution is obliged to process the product before its disposal, see Care and handling.*

*Always adhere to national regulations when disposing of or recycling the product or its components!*

Products carrying this symbol are subject to separate collection of electric and electronic devices. Within the European Union, disposal is carried out without fees by the manufacturer.

If you have any questions concerning the disposal of the product, please contact your national B. Braun/Aesculap agency, see Technical Service.

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Aesculap AG