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Read the instructions for use carefully before working with the unit.

Only qualified personnel may use the saw together with accessories.

Keep the instructions for use safely for later use. They must be available to the theatre staff at all times.

We reserve the right to make changes to the models and technical specifications without prior notice.

Combinations of accessories other than those stated in the instructions for use may be used only if they do not impair the features and safety requirements for the intended applications.

Only Aesculap products may be combined!

1. Overview

- 1 ACCULAN sagittal saw GA 633
- 2 Closure GA 624
- 3 Closure release
- 4 Battery GA 616
- 5 Loading side
- 6 Funnel GA 625
- 7 Saw blade
- 8 Trigger for controlling the oscillation frequency
- 9 Sliding catch on knob/pushbutton
- 10 Knob/pushbutton for fixing saw blade
- 11 Saw head
- 12 Sliding catch on tongue
- 13 Tongue for locking saw head
- 14 Changing funnel

GA 633 delivery

ACCULAN sagittal saw	GA 633
Changing funnel	GA 614
Closure	GA 624
Funnel	GA 625
Battery pack	GA 616
Instructions for use	TA 010 010
Storage tray (for GA 633)	GB 509 R

Explanation of symbols

	Attention, see instructions for use
	CE marking according to directive 93/42/EEC
	Opening/Closing
	Classification type B
	Control of oscillation frequency: • Trigger not pressed (0): machine not operating (0) • The oscillation frequency increases as the trigger is depressed further. • Trigger completely depressed (I): machine running at maximum speed (I)
	Protect from moisture!
	Marking of electric and electronic devices according to directive 2002/96/EC (WEEE), see Disposal.

Table of Contents

1. Overview	11
2. Intended use	12
3. Mode of operation	12
4. Danger information and warnings	12
5. Operation	12
6. Setup	13
6.1 Battery insertion (Fig. 1-3)	13
6.2 Unsterile battery removal (Fig. 4)	14
6.3 Sterile battery removal (Fig. 5)	14
7. Operation	15
7.1 Fixing the saw blade	15
7.2 Releasing the saw blades	16
7.3 Turning the saw head	16
8. Functional test	16
9. Care and handling	17
9.1 Manual preparation	17
9.2 Mechanical preparation	17
10. Notes on use of battery GA 616	18
11. Recognizing and remedying errors	18
12. Maintenance	19
13. Technical Service	19
14. Technical data	19
15. Accessories/spare parts	19
16. Disposal	19

2. Intended use

The ACCULAN sagittal saw GA 633 is used for cutting bones in orthopaedics/traumatology.

3. Mode of operation

In the ACCULAN sagittal saw GA 633, the electric motor is supplied with voltage by battery. The oscillation frequency is electronically controlled. With the trigger **8**, this oscillation frequency can be continuously controlled up to 15,000 ¹/min.

Note

*The continuous regulation of the oscillation frequency necessitates electronic components. The conditions (steam, temperature, etc.) prevailing during the sterilization impose a considerable load on these components. To avoid this load, the complete electronics of the battery saw is integrated in the battery **4**, which may not be sterilized. The battery **4** too should therefore be handled with the care usual for a surgical instrument.*

Bones are cut as a result of the to-and-fro movement of the saw blade. The small stroke prevents damage to the soft-tissue.

No spanners are required for tensioning the saw blades.

You can turn the saw head of the sagittal saw to four different positions.

4. Danger information and warnings

WARNING

Inadvertent activation of the drive can lead to injury from the saw blade!

- Only replace the saw blades once the drive has been secured against inadvertent activation!
- Turn the catch on the trigger **8** to the OFF position. Also see under Section 7.

CAUTION

Blunt saw blades and poorly maintained or unmaintained machines can quickly heat up and may cause burns on skin and tissue!

CAUTION

- ACCULAN sagittal saw GA 633 may only be steam-sterilized!

CAUTION

The batteries GA 616 may on no account be sterilized since heating to above 65 °C causes irreparable damage.

- For notes on the use of battery GA 616, see under Section 10.

CAUTION

The ACCULAN sagittal saw and the accessories are delivered in the unsterile state.

- Before use the saw and the accessories must be prepared as detailed in Section 9.

Note

Only Aesculap products may be combined!

5. Operation

The following components are required for putting the ACCULAN sagittal saw GA 633 into operation:

- ACCULAN charger GA 627 M
- ACCULAN battery pack GA 616
- ACCULAN sagittal saw GA 633

- Before putting the saw into operation for the first time, ensure that the required accessories are available.
- Check the individual parts for visible damage.
- Use only satisfactory machines and accessories.

CAUTION

The ACCULAN sagittal saw and the accessories are delivered in the unsterile state.

- Before being used for the first time, the battery saw GA 633 and the accessories are prepared as described under Section 9.

Note

Consult the special instructions for use from Aesculap for putting the ACCULAN charger into operation.

Please also see the separate notes on the handling and care of the battery GA 616 under Section 10.

The ACCULAN sagittal saw GA 633 can be operated with saw blades with the appropriate connection (see Section 14).

6. Setup

6.1 Battery insertion (Fig. 1–3)

- Products that are not marked as sterile must be sterilized according to the respective instructions prior to use.

CAUTION

- The batteries GA 616 may on no account be sterilized since heating to above 65 °C causes irreparable damage.

The following procedure should be followed for the sterile preparation of the battery saw:

- The funnel **6** (sterile) is mounted on the machine held upside down (battery shaft facing upwards).


Fig. 1

- The battery **4** (unsterile) is then introduced into the charging shaft by a second person (the battery cannot rotate), and the funnel **6** is removed and is placed in the unsterile room.


Fig. 2

- The closure **2** (sterile) is then mounted and is securely snapped in. The closure ensures complete tightness so that there is no danger of contamination of the operating site by the unsterile battery.


Fig. 3

6.2 Unsterile battery removal (Fig. 4)

After the end of the operation, the battery must be removed prior to cleaning and disinfection and especially before sterilization.

- For this purpose, you must simultaneously press the two release buttons **3** on the closure **2** to enable you to remove it (the closure cannot be opened if only one release button is activated). This ensures that inadvertent opening by pressing only one side is prevented.
- Prior to battery removal, make sure that the trigger for controlling the oscillation frequency **8** is in the front position (rest position).

CAUTION

Tapping on hard objects damages the machine and also the battery.

- By tapping on the palm of the hand or on a soft surface, the battery can be removed from the battery shaft.

Fig. 4

6.3 Sterile battery removal (Fig. 5)

The changing funnel GA 614 is used for intraoperative battery change while safely maintaining the sterile conditions. Procedure:

- Remove the cover **2** of the machine.
- Prior to battery removal, make sure that the trigger for controlling the oscillation frequency **8** is in the front position (rest position).
- Mount the sterile changing funnel **14** on the machine.
- Gentle shaking/tapping causes the battery to slide gently into the changing funnel.
- The changing funnel including the flat battery is handed to the unsterile person.
- The ACCULAN machine is equipped again with a fully charged battery as described in Section 6.1.

Fig. 5

7. Operation

Prevention of inadvertent activation

To prevent inadvertent activation (e.g. when changing the saw blade), the machine is equipped with a catch on the trigger for controlling the oscillation frequency.

Explanation of symbols

ON: The machine is ready for operation.

OFF: The trigger **8** is blocked. The machine cannot be put into operation.

- a) By turning the trigger catch to the OFF position, you can block the trigger for controlling the oscillation frequency. Inadvertent operation of the machine is thus prevented.
- b) By turning the trigger catch to the ON position, the trigger for controlling the oscillation frequency is released. The machine is now ready for operation.


Fig. 6

- When the trigger **8** is pressed, the machine operates with a controllable oscillation frequency of 0-15,000 $1/\text{min}$.

Note

Owing to the electronics used, a quiet whistling sound is heard shortly before the soft start-up of the sagittal saw; however, this sound increasingly diminishes immediately the saw starts operating again.

- Begin with a moderate oscillation frequency until the saw blade is guided in the cut and then increase the oscillation frequency during the sawing process.
- The best cutting performance is achieved at high drive speed in combination with a slight to-and-fro movement of the sagittal saw; this results in better chip removal.
- Take care not to bend or twist the saw blade during the sawing process. Otherwise, there is a danger that the saw blade will break.

7.1 Fixing the saw blade

WARNING

Inadvertent activation of the saw can cause injuries from the saw blade!

- Ensure that the trigger is blocked before changing the saw blade!
- Turn the catch on the trigger **8** to the OFF position.

Preparation

- Push the sliding catch **9** to the left.
- Turn the knob/pushbutton **10** several revolutions counterclockwise as far as it will go.

Fig. 7

Insert saw blade

- Press the knob/pushbutton **10** towards the saw head.
- Insert the saw blade **7** into the pins of the saw blade coupling.


Fig. 8

Fix saw blade

- Press the sliding catch **9** to the right.
- Turn the rotary knob/pushbutton **10** clockwise to the stop.
- Tighten the rotary knob/pushbutton **10**.
The saw blade coupling is now secured to prevent inadvertent opening.


Fig. 9

Functional test

- Check whether the saw blade fits firmly in the coupling by bending the saw blade up and down with gentle pressure. The saw blade coupling must not open.
- The sliding catch **9** must have been pushed fully to the right.

7.2 Releasing the saw blades

- Slide the sliding catch **9** to the left and fully undo the rotary/push-button **10** by turning anti-clockwise.
- Press the knob/pushbutton towards the saw head.
- Remove the saw blade **7** from the saw blade coupling.

Fig. 10

7.3 Turning the saw head

WARNING

Inadvertent activation of the saw can cause injuries from the saw blade.

- Make sure that the trigger is blocked before turning the saw head!
- Turn the catch on the trigger **8** to the OFF position.

Note

The saw may be activated only when the saw head is locked to prevent turning. Working with inadvertently released saw head may cause damage to the sagittal saw.

You can turn the saw head **11** to four different positions:


Fig. 11

- The saw head **11** cannot be turned through 360° but in each case only up to a stop.
- The saw head must snap in perceptibly in each of the four positions.
- Before sawing, the saw head must be locked in the respective position.

Preparation

- To turn the saw head **11**, push the sliding catch **12** to the left.
- Turn the tongue **13** half a revolution counterclockwise.

Fig. 12

Turning and locking

- The saw head **11** can now be turned through 90° to the desired position.
- To lock the saw head, push the sliding catch **12** to the right and turn the tongue **13** half a revolution clockwise as far as it will go.


Fig. 13

8. Functional test

Each time before use, you must subject the ACCULAN sagittal saw with the inserted saw blade to a test run.

- To do this, first check whether the saw blade is firmly secured in the coupling (see Section 7).
- Operate the sagittal saw briefly at maximum speed. Check that it is functioning properly and reliably. The sagittal saw must not heat up to an unusual extent nor make unusual noises.
- Each time before use, check the saw blades in the unclamped state for damage and wear.

Note

Each intraoperative battery change has to be followed by a functional test!

CAUTION

Blunt saw blades and poorly maintained or unmaintained machines can heat up very quickly and may cause burns on skin and tissue!

- Only use fault-free saw blades.
- Replace blunt or damaged saw blades.

9. Care and handling

9.1 Manual preparation

9.1.1 Cleaning and disinfecting

- The ACCULAN sagittal saw can be cleaned and disinfected with all commercially available cleaning and disinfecting solutions that are suitable for cleaning anodized aluminum surfaces and that may be used in conformance with national and state laws. The manufacturer's instructions with regard to dilution ratios, action time, etc. must be followed.
- Clean the ACCULAN sagittal saw with a soft brush or with a cloth that has been moistened in disinfecting solution.

CAUTION

- **Do not lay, immerse or heat the ACCULAN sagittal saw in any fluid.** Laying or immersing the saw in fluid will damage its mechanical components. Should any fluid penetrate the unit, allow it to run out immediately.

9.1.2 Maintenance

The ACCULAN sagittal saw requires no maintenance.

9.1.3 Sterilization

- In order to achieve safe and reliable sterilization of the saw, the closure **2** must first be removed.
- Sterilize the ACCULAN sagittal saw **1**, closure **2** funnel **6** and changing funnel **14** in a steam sterilizer only.
- Sterilize with steam, taking note of the following:
Carry out steam sterilization via an approved steam sterilization procedure (e.g. sterilizer in conformance with EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993 and validated in conformance with EN 554/ISO 13683). Carry out prevacuum sterilization for a minimum of 5 minutes at 27 °F/134 °C and at 2 bar pressure.
- As needed, lay the components on a textile surface or wrap them in cloths.
- Position the ACCULAN sagittal saw in the sterilizer in such a way that any condensation that forms in the battery shaft will run out. The battery shaft must be absolutely dry after sterilization.

CAUTION

- **Do not under any circumstances sterilize the rechargeable batteries, since heating them to a temperature above 65 °C will cause irreparable damage to them.**
- For information on the use of the rechargeable battery GA 616, see Section 10.

9.2 Mechanical preparation

9.2.1 Cleaning and disinfecting

CAUTION

Improper mechanical cleaning will damage the ACCULAN sagittal saw!

- Clean the ACCULAN sagittal saw in a suitable storage tray only.
- Use a suitable cleaning and disinfecting cycle for mechanical cleaning.
- Do not use alkaline cleaning agents.
- Perform mechanical cleaning immediately after the operation.
- Place storage tray GB 509 R in a suitable wire basket (e.g. JF 214 R).
- Position the ACCULAN sagittal saw properly in storage tray GB 509 R immediately after the operation.
- To avoid abrupt changes in temperature, use a single-chamber washer disinfectant.
- Select a cleaning and disinfecting cycle according to the following criteria:
 - Neutral cleaner
 - Cleaning temperature ≤ 50 °C
 - Final rinse with simultaneous thermal disinfection in demineralized water.
- Example: Miele Vario TD process.
- Follow the instructions for use of the Aesculap ECCOS storage system TA 009 721.

9.2.2 Maintenance

The ACCULAN sagittal saw requires no maintenance.

9.2.3 Sterilization

Sterilize the ACCULAN sagittal saw immediately after it has been cleaned mechanically.

- The saw can remain in the storage tray during sterilization.
- Sterilize with steam, taking note of the following:
Carry out steam sterilization via an approved steam sterilization procedure (e.g. sterilizer in conformance with EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993 and validated in conformance with EN 554/ISO 13683). Carry out prevacuum sterilization for a minimum of 5 minutes at 274 °F/134 °C and at 2 bar pressure.

CAUTION

- **Do not under any circumstances sterilize the rechargeable batteries, since heating them to a temperature above 65 °C will cause irreparable damage to them.**
- For information on the use of the rechargeable battery GA 616, see Section 10.

10. Notes on use of battery GA 616

The GA 616 battery is a high-performance NiMH battery with integrated electronics and power supply. Thanks to its unique design and computer-controlled charger, the battery has an unusually long service life.

IMPORTANT

Heating the rechargeable battery to temperatures exceeding 65 °C will reduce the battery's capacity. Subjecting the battery to such temperatures for an extended period will ruin it.

- Do not heat the battery to temperatures exceeding 65 °C.

IMPORTANT

If the poles of the battery (either top side or down side) are connected by metal or liquid, high current will flow (comparable with a car battery). The battery may be heated up extremely and be destroyed and it may even damage surrounding objects (danger of fire).

- Avoid a short-circuit.

Therefore we recommend to treat the battery as follows:

- Store the battery at room temperature.

- Do not sterilize or disinfect the battery by any means. As the battery is completely sealed inside the housing of the machine the contamination of the OR-field by the battery is not possible. Thus, disinfection or sterilization is not necessary.
- Always keep the battery dry.
- Never immerse the battery in fluids.
- Never place the battery on wet surfaces.
- Always make sure that the sockets of the charging unit are dry inside before you put the battery in. If this is not the case, dry the sockets with a cloth.
- Do not place or run the charging unit and the battery near to places where you handle fluids (instrument cleaning, infusions, etc.).

If you take these measures, your ACCULAN system is sure to serve you well. If you have any questions concerning the ACCULAN system, contact your medical product adviser.

Protect from moisture!

11. Recognizing and remedying errors

Fault	Cause	Fault detection	Remedy
Saw blade not moving	Defective gearing	Motor makes audible sound when running	Repair by Aesculap service department
Inadequate cutting performance of saw blades	Blunt saw blade	Worn saw blade teeth	Insert new saw blade
	Performance or speed of drive insufficient	NiMH battery is too weak NiMH battery is defective or flat	Charge battery Replace battery
	Inadequate chip removal	Hot saw blade	To-and-fro movement during sawing
Saw blade cannot be coupled	Dirt	Dirt on the saw blade and in the coupling	Clean coupling
	Deformation	Deformed coupling parts or saw blades	Repair by Aesculap service department
	Wrong saw blade	Saw blade not from current range	Insert correct saw blade (see Section 14)
	Rotary knob/pushbutton not opened and/or not depressed	Rotary knob/pushbutton not turned fully anti-clockwise	Couple saw blade according to Section 7.1
Saw blade cannot be uncoupled	Rotary knob/pushbutton not opened	Rotary knob/pushbutton not turned fully anti-clockwise	Uncouple saw blade according to Section 7.2

12. Maintenance

To ensure reliable operation, Aesculap recommends that maintenance be performed after 300 reprocessing cycles, at least once per year.

- For services to this end, please contact an Aesculap/B. Braun representative in your country, see Technical Service.

Note

The battery compartment of the machine is closed germ-tight and water-tight with the closure.

- You should therefore regularly check the silicone seal in the closure and the sealing surface on the handle of the saw. A damaged seal and/or sealing surface can result in contamination with germs from the battery shaft.
- Replace blunt or damaged saw blades.

CAUTION

Blunt saw blades and poorly maintained or unmaintained machines can heat up very quickly and may cause burns on skin and tissue!

13. Technical Service

For service, maintenance or repairs contact your Aesculap/B. Braun representative.

Carrying out modifications on medical equipment renders the guarantee, and any licenses, null and void.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-2700
Fax: +49 7461 16-2887
E-mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

14. Technical data

Type of unit	GA 633
DC voltage	9.6 V
Rated capacity of the battery	1.6 Ah
Rated current	8 A
Max. mechanical power	120 W
Oscillation frequency	0–15,000 ¹ /min
Dimensions	170 mm x 200 mm x 45 mm
Weight (with battery)	approx. 1,800 g
Type (according to IEC 601)	B
Electromagnetic compatibility	IEC 601-1-2
Standards applied	IEC 601-1

15. Accessories/spare parts

Charger	GA 627 M
Battery	GA 616
Closure	GA 624
Sterile funnel	GA 625
Changing funnel	GA 614
Saw blade	e.g. GC 677 R
Storage tray (for GA 633)	GB 509 R
Storage tray (for GA 624)	GB 504 R
Storage tray (for GA 614)	GB 506 R

16. Disposal

Note

The user institution is obliged to process the product before its disposal, see Care and handling.

Always adhere to national regulations when disposing of or recycling the product or its components!

Products carrying this symbol are subject to separate collection of electric and electronic devices. Within the European Union, disposal is carried out without fees by the manufacturer.

If you have any questions concerning the disposal of the product, please contact your national B. Braun/Aesculap agency, see Technical Service.

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Aesculap AG