# DK50 PLUS DK50 2V

- EN User manual
- DE Benutzerhandbuch
- FR Manuel de l'utilisateur
- RU Руководство пользователя
- PL Instrukcja obsługi
- SK Návod na použitie
- CS Návod k použití

COMPRESSOR  
KOMPRESSOR  
COMPRESSEUR  
КОМПРЕССОР  
SPREŻARKA  
KOMPRESOR

DK50 PLUS  
DK50 2V

**EKOM spol. s r. o.**  
Priemyselná 5031/18  
SK-921 01 Piešťany  
Slovak Republic  
tel.: +421 33 7967255  
fax: +421 33 7967223

[www.ekom.sk](http://www.ekom.sk)  
email: [ekom@ekom.sk](mailto:ekom@ekom.sk)

DATE OF LAST REVISION  
DATUM DER LETZTEN ÜBERARBEITUNG  
DATE DE LA DERNIÈRE RÉVISION  
ДАТА ПОСЛЕДНЕЙ РЕДАКЦИИ  
DATA OSTATNIEJ AKTUALIZACJI  
DÁTUM POSLEDNEJ REVÍZIE  
DATUM POSLEDNÍ REVIZE

12/2018

NP-DK50 PLUS, 2V-5\_12-  
2018-MD  
112000359-000

EN	<u>CONTENTS..... 1</u>
DE	<u>INHALT..... 39</u>
FR	<u>TABLES DES MATIÈRES..... 78</u>
RU	<u>ОГЛАВЛЕНИЕ.....117</u>
PL	<u>SPIS TREŚCI.....156</u>
SK	<u>OBSAH.....195</u>
CS	<u>OBSAH.....234</u>

## CONTENTS

<b>IMPORTANT INFORMATION .....</b>	<b>2</b>
1. CONFORMITY WITH THE REQUIREMENTS OF EUROPEAN UNION DIRECTIVES .....	2
2. INTENDED USE .....	2
3. CONTRAINDICATIONS AND SIDE-EFFECTS.....	2
4. WARNINGS AND SYMBOLS .....	2
5. WARNINGS .....	3
6. STORAGE AND TRANSPORT .....	5
7. TECHNICAL DATA .....	6
8. PRODUCT DESCRIPTION .....	12
9. PRODUCT FUNCTIONALITY .....	14
10.PNEUMATIC DIAGRAM .....	18
<b>INSTALLATION .....</b>	<b>19</b>
11.CONDITIONS FOR USE .....	19
12.PLACEMENT OF THE COMPRESSOR .....	20
13. PNEUMATIC CONNECTIONS .....	23
14. ELECTRICAL CONNECTIONS .....	24
15. CONNECTION DIAGRAM .....	26
<b>OPERATION .....</b>	<b>29</b>
16.COMMISSIONING .....	29
17.SWITCHING THE COMPRESSOR ON .....	30
18.COMPRESSOR SHUT-DOWN .....	30
<b>MAINTENANCE .....</b>	<b>31</b>
19.DEVICE MAINTENANCE .....	31
<b>TROUBLESHOOTING .....</b>	<b>37</b>
20.REPAIR SERVICE .....	38
21.STORAGE .....	38
22.DISPOSAL OF THE DEVICE .....	38
<b>ANNEX .....</b>	<b>273</b>
23.INSTALLATION RECORD.....	273

## **IMPORTANT INFORMATION**

### **1. CONFORMITY WITH THE REQUIREMENTS OF EUROPEAN UNION DIRECTIVES**

This product conforms to the requirements of the Medical Device Directive (93/42/EEC) and the Machinery Directive (2006/42/EC), and is safe for the intended use if all safety instructions are followed.

### **2. INTENDED USE**

The compressor is used as a source of clean, oil-free compressed air to power dental unit, instruments, and devices in laboratories where the parameters and properties of the compressed air are suitable for the specific application.

**Compressed air supplied by the compressor is unsuitable for use with artificial lung ventilation devices without further filtration**

Any other use of the product beyond its intended use is considered an incorrect use. The manufacturer is not liable for any damages or injuries resulting from the incorrect use.

### **3. CONTRAINDICATIONS AND SIDE-EFFECTS**

There are no contraindications or side-effects known.

### **4. WARNINGS AND SYMBOLS**

The following symbols are used in the user manual, device and its packaging to denote important details and information:

General warnings

Warning

Danger, electric shock hazard

Read the user manual!

CE-marking

Compressor is controlled automatically and may start without warning

Caution! Hot surface

Ground connection

Terminal for ground connection

Alternating current

Handling mark on package – FRAGILE

Handling mark on package – THIS SIDE UP

Handling mark on package – KEEP DRY

Handling mark on package – TEMPERATURE LIMITATIONS

Handling mark on package – LIMITED STACKING

Mark on package – RECYCLABLE MATERIAL

Manufacturer

## 5. WARNINGS

The product is designed and manufactured to be safe for the user and the surrounding environment when used in the defined manner. Keep the following warnings in mind.

### 5.1. General warnings

**READ THE USER MANUAL CAREFULLY BEFORE USING THE DEVICE AND KEEP IT FOR FUTURE USE!**

- The user manual supports the correct installation, operation and maintenance of the product. Careful review of this manual will provide the information necessary for the correct operation of the product for its intended use.
- Original packaging should be retained for the return of the device. Only the original packaging ensures protection of the device during transport. If it is necessary to return the product during the warranty period, the manufacturer is not liable for damages caused by improper packaging.
- This warranty does not cover damages originating from the use of accessories other than those specified or recommended by the manufacturer.
- The manufacturer only guarantees the safety, reliability and function of the device if:
  - installation, new settings, changes, expansion, and repairs are performed by the manufacturer or an organisation authorised by the manufacturer.
  - the product is used pursuant to the user manual.

- The user manual corresponds to the configuration of the product and its compliance with the applicable safety and technical standards at the time of printing. The manufacturer reserves all rights for the protection of its configuration, methods and names.
- Translation of the user manual is performed in accordance with the best available knowledge. The Slovak version is to be used in the event of any uncertainties.

## 5.2. General safety warnings

The manufacturer designed and manufactured the product to minimise all risks when used correctly for the intended use. The manufacturer considers it its obligation to lay down the following general safety precautions.

- Use and operation of the product must comply with all laws and local regulations valid in the place of use. The operator and user are responsible for following all the appropriate regulations in the interests of performing work safely.
- Only the use of original parts guarantees the safety of operating personnel and the flawless operation of the product itself. Only the accessories and parts mentioned in the technical documentation or expressly approved by the manufacturer should be used.
- The operator must ensure that the device is functioning correctly and safely before every use.
- The user must be familiar with the operation of the device.
- Do not use the product in environments with a risk of explosion.
- The user must inform the supplier immediately if any problem directly related to the operation of the device occurs.

## 5.3. Safety warnings on protection from electric current

- The device must only be connected to a properly installed, earthed socket.
- Before the product is plugged in, ensure that the mains voltage and frequency stated on the product are in compliance with the values of the mains.
- Prior to putting the device into operation check for any damage to the connected pneumatic lines and electrical wiring. Replace damaged pneumatic lines and electrical wirings immediately.
- Immediately disconnect the product from the mains (remove power cord from the socket) in hazardous situations or in the case of a technical malfunction.
- During all repairs and maintenance, ensure that:
  - the mains plug is removed from the power socket
  - pressure is vented from the air tank and pipes
- The product shall only be installed by a qualified technician.

## 6. STORAGE AND TRANSPORT

The compressor is shipped from the factory in a transport packaging. This protects the device from damage during transport.

**The original compressor packaging must be used for transport whenever possible. Transport the compressor in an upright position, always secured with transport fixation.**

Protect the compressor from moisture, dirt and extreme temperatures during transport and storage. Store the compressor in its original packaging in a warm, dry, and dust-free area. Do not store near any chemical substances.

Keep the packaging material, if possible. If not, please dispose of the packaging material in an environmentally-friendly way. Packaging cardboard can be recycled with old paper.

**The compressor may only be transported when all air pressure has been vented. Before moving or transporting the compressor, release all air pressure from the tank and hoses and drain condensate from the air tank.**

**Storing or shipping the equipment in any conditions other than those specified below is prohibited.**

### Ambient conditions for storage and transport

Products may only be stored and transported in vehicles that are free of any traces of volatile chemicals under the following conditions:

Temperature:  $-25^{\circ}\text{C}$  to  $+55^{\circ}\text{C}$ , 24 h at up to  $+70^{\circ}\text{C}$

Relative humidity: 10% to 90% (non-condensing)

## 7. TECHNICAL DATA

Compressors are designed for dry and ventilated indoor environments with the following conditions:

Temperature : +5°C to +40°C,  
 Max. relative humidity.: 70%,  
 Max. absolute humidity. 15 g/m<sup>3</sup>

Tab.1

5 – 7 bar		DK50 PLUS		DK50 PLUS S		DK50 PLUS/M		DK50 PLUS S/M	
Nominal voltage / frequency (*)	V / Hz	230 / 50/60	115 / 60	230 / 50/60	115 / 60	230 / 50/60	115 / 60	230 / 50/60	115 / 60
Performance at 5 bar	Lit.min <sup>-1</sup>	75/85	85	75/85	85	58/68	68	58/68	68
Working pressure (**)	bar	5.0 – 7.0		5.0 – 7.0		5.0 – 7.0		5.0 – 7.0	
Performance with KJF-1 at 5 bar	Lit.min <sup>-1</sup>	75/85	85	75/85	85	-		-	
Max. current	A	3.8/4.8	8.9	3.9/4.9	9	3.9/4.9	9	4/5	9.1
Motor power	kW	0.55		0.55		0.55		0.55	
Air tank volume	Lit.	25		25		25		25	
Air quality – filtration	µm	-		-		0.3		0.3	
Maximum operating pressure of safety valve	bar	8.0		8.0		8.0		8.0	
Noise level at 5 bar	LpfA [dB]	≤ 66/68	≤ 68	≤ 47/50	≤ 50	≤ 66/68	≤ 68	≤ 47/50	≤ 50
Duty cycle		100%		100%		100%		100%	
PDP - drying performance at 7 bar		-		-		≤ +3°C		≤ +3°C	
Time to fill air tank from 0 to 6 bar	s	123/105	105	123/105	105	157/128	128	157/128	128
Dimensions (net) w x d x h	mm	460x500x708		560x675x877		460x565x708		560x675x877	
Net weight (****)	kg	48 (***)		87 (***)		53		92	
Classification under EN 60601-1		Class I.							



Notes:

(\*) When ordering, state the version of the compressor

(\*\*) For other range of pressure:consult with the supplier

(\*\*\*) Weight of compressors with a KJF-1 unit \_add 3 kg

(\*\*\*\*) Information about the weight is for informative purposes only and applies to the product without any additional accessories.

Tab. 2

6 – 8 bar		DK50 PLUS		DK50 PLUS S		DK50 PLUS/M		DK50 PLUS S/M	
Nominal voltage / frequency (*)	V / Hz	230 / 50/60	115 / 60	230 / 50/60	115 / 60	230 / 50/60	115 / 60	230 / 50/60	115 / 60
Performance at 6 bar	Lit.min <sup>-1</sup>	70/80	80	70/80	80	60/70	70	60/70	70
Working pressure (**)	bar	6.0 – 8.0		6.0 – 8.0		6.0 – 8.0		6.0 – 8.0	
Performance with KJF-1 at 6 bar	Lit.min <sup>-1</sup>	70/80	80	70/80	80	-		-	
Max. current	A	3.9/4.9	9	4/5	9.1	4/5	9.1	4.1/5.1	9.2
Motor power	kW	0.55		0.55		0.55		0.55	
Air tank volume	Lit.	25		25		25		25	
Air quality – filtration	µm	-		-		0.3		0.3	
Maximum operating pressure of safety valve	bar	12.0		12.0		12.0		12.0	
Noise level at 5 bar	LpFA [dB]	≤ 66/68	≤ 68	≤ 47/50	≤ 50	≤ 66/68	≤ 68	≤ 47/50	≤ 50
Duty cycle		100%		100%		100%		100%	
PDP drying performance at 7 bar		-		-		≤ +3°C		≤ +3°C	
Time to fill air tank from 0 to 7 bar	s	149/126	126	149/126	126	198/159	159	198/159	159
Dimensions (net) w x d x h	mm	460x500x708		560x675x877		460x565x708		560x675x877	
Net weight (****)	kg	48 (***)		87 (***)		53		92	
Classification under EN 60601-1		Class I.							



## Notes:

(\*) When ordering, state the version of the compressor

(\*\*) For other range of pressure consult with the supplier

(\*\*\*) Weight of compressors with a KJF-1 unit \_add 3 kg

(\*\*\*\*) Information about the weight is for informative purposes only and applies to the product without any additional accessories.

Tab. 3

8 – 10 bar		DK50 PLUS		DK50 PLUS S		DK50 PLUS/M		DK50 PLUS S/M	
Nominal voltage / frequency (*)	V / Hz	230 / 50/60	115 / 60	230 / 50/60	115 / 60	230 / 50/60	115 / 60	230 / 50/60	115 / 60
Performance at 8 bar	Lit.min <sup>-1</sup>	60/70	70	60/70	70	50/60	60	50/60	60
Working pressure (**)	bar	8.0 – 10.0		8.0 – 10.0		8.0 – 10.0		8.0 – 10.0	
Performance with KJF-1 at 8 bar	Lit.min <sup>-1</sup>	60/70	70	60/70	70	-		-	
Max. current	A	4.1/5.1	9.2	4.2/5.2	9.3	4.2/5.2	9.3	4.3/5.3	9.4
Motor power	kW	0.55		0.55		0.55		0.55	
Air tank volume	Lit.	25		25		25		25	
Air quality – filtration	µm	-		-		0.3		0.3	
Maximum operating pressure of safety valve	bar	12.0		12.0		12.0		12.0	
Noise level at 5 bar	Lp <sub>f</sub> A [dB]	≤ 66/68	≤ 68	≤ 47/50	≤ 50	≤ 66/68	≤ 68	≤ 47/50	≤ 50
Duty cycle		100%		100%		100%		100%	
PDP drying performance at 7 bar		-		-		≤ +3°C		≤ +3°C	
Time to fill air tank from 0 to 9 bar	s	208/176	176	208/176	176	304/239	239	304/239	239
Dimensions (net) w x d x h	mm	460x500x708		560x675x877		460x565x708		560x675x877	
Net weight (****)	kg	48 (***)		87 (***)		53		92	
Classification under EN 60601-1		Class I.							



## Notes:

(\*) When ordering, state the version of the compressor

(\*\*) For other range of pressure consult with the supplier

(\*\*\*) Weight of compressors with a KJF-1 unit \_add 3 kg

(\*\*\*\*) Information about the weight is for informative purposes only and applies to the product without any additional accessories.

Tab.4

5 – 7 bar		DK50 2V		DK50 2VS		DK50 2V/M		DK50 2VS/M	
Nominal voltage / frequency (*)	V / Hz	230/50 3x400/50	230/60 115/60	230/50 3x400/50	230/60 115/60	230/50 3x400/50	230/60 115/60	230/50 3x400/50	230/60 115/60
Performance at 5 bar	Lit.min <sup>-1</sup>	140		140		108		108	
Working pressure (**)	bar	5.0 – 7.0		5.0 – 7.0		5.0 – 7.0		5.0 – 7.0	
Performance with KJF at 5 bar	Lit.min <sup>-1</sup>	140		140		-		-	
Max. current	A	8.3 3.6	9.0 18.2	8.5 3.8	9.2 18.4	8.4 3.7	9.1 18.3	8.6 3.9	9.3 18.5
Motor power	kW	1.2 1.2	1.2 1.1	1.2 1.2	1.2 1.1	1.2 1.2	1.2 1.1	1.2 1.2	1.2 1.1
Air tank volume	Lit.	25		25		25		25	
Air quality – filtration	µm	-		-		0,3		0,3	
Maximum operating pressure of safety valve	bar	8.0		8.0		8.0		8.0	
Noise level at 5 bar	LpA [dB]	≤ 67		≤ 51.5		≤ 69.0		≤ 55.0	
Duty cycle		100%		100%		100%		100%	
PDP - drying performance at 7 bar		-		-		≤ +3°C		≤ +3°C	
Time to fill air tank from 0 to 6 bar	s	52		52		65		65	
Dimensions (net) w x d x h	mm	460x500x709		560x675x877		460x565x709		560x675x877	
Net weight (****)	kg	55(***) 53(***)	55(***)	93(***) 91(***)	93(***)	60 58	60	98 96	98
Classification under EN 60601-1		Class I.							



## Notes:

(\*) When ordering, state the version of the compressor

(\*\*) For other range of pressure:consult with the supplier

(\*\*\*) Weight of compressors with a KJF-1(KJF-2) unit \_add 3(3.5) kg

(\*\*\*\*) Information about the weight is for informative purposes only and applies to the product without any additional accessories

Tab. 5

6 – 8 bar		DK50 2V		DK50 2VS		DK50 2V/M		DK50 2VS/M	
Nominal voltage / frequency (*)	V / Hz	230/50 3x400/50	230/60 115/60	230/50 3x400/50	230/60 115/60	230/50 3x400/50	230/60 115/60	230/50 3x400/50	230/60 115/60
Performance at 6 bar	Lit.min <sup>-1</sup>	135		135		104		104	
Working pressure (**)	bar	6.0 – 8.0		6.0 – 8.0		6.0 – 8.0		6.0 – 8.0	
Performance with KJF at 6 bar	Lit.min <sup>-1</sup>	135		135		-		-	
Max. current	A	8.6 3.8	9.2 19.9	8.8 4.0	9.4 20.1	8.7 3.9	9.3 20.0	8.9 4.1	9.5 20.2
Motor power	kW	1.2 1.2	1.2 1.1	1.2 1.2	1.2 1.1	1.2 1.2	1.2 1.1	1.2 1.2	1.2 1.1
Air tank volume	Lit.	25		25		25		25	
Air quality – filtration	µm	-		-		0.3		0.3	
Maximum operating pressure of safety valve	bar	12.0		12.0		12.0		12.0	
Noise level at 5 bar	LpA [dB]	≤ 67		≤ 51.5		≤ 69.0		≤ 55.0	
Duty cycle		100%		100%		100%		100%	
PDP drying performance at 7 bar		-		-		≤ +3°C		≤ +3°C	
Time to fill air tank from 0 to 7 bar	s	61		61		78		78	
Dimensions (net) w x d x h	mm	460x500x720		560x675x877		460x565x720		560x675x877	
Net weight (****)	kg	55 <sup>(***)</sup> 53 <sup>(***)</sup>	55 <sup>(***)</sup>	93 <sup>(***)</sup> 91 <sup>(***)</sup>	93 <sup>(***)</sup>	60 58	60	98 96	98
Classification under EN 60601-1		Class I.							



## Notes:

(\*) When ordering, state the version of the compressor

(\*\*) For other range of pressure consult with the supplier

(\*\*\*) Weight of compressors with a KJF-1(KJF-2) unit \_add 3(3.5) kg

(\*\*\*\*) Information about the weight is for informative purposes only and applies to the product without any additional accessories.

Tab. 6

8 – 10 bar		DK50 2V		DK50 2VS		DK50 2V/M		DK50 2VS/M	
Nominal voltage / frequency (*)	V / Hz	230/50 3x400/50	230/60 115/60	230/50 3x400/50	230/60 115/60	230/50 3x400/50	230/60 115/60	230/50 3x400/50	230/60 115/60
Performance at 8 bar	Lit.min <sup>-1</sup>	110		110		87		87	
Working pressure (**)	bar	8.0 – 10.0		8.0 – 10.0		8.0 – 10.0		8.0 – 10.0	
Performance with KJF at 8 bar	Lit.min <sup>-1</sup>	110		110		-		-	
Max. current	A	8.4 4.6	9.3 18.6	8.6 4.8	9.5 18.8	8.5 4.7	9.4 18.7	8.7 4.9	9.6 18.9
Motor power	kW	1.2 1.2	1.2 1.1	1.2 1.2	1.2 1.1	1.2 1.2	1.2 1.1	1.2 1.2	1.2 1.1
Air tank volume	Lit.	25		25		25		25	
Air quality – filtration	µm	-		-		0.3		0.3	
Maximum operating pressure of safety valve	bar	12,0		12,0		12,0		12,0	
Noise level at 5 bar	LpA [dB]	≤ 67.0		≤ 51.5		≤ 69.0		≤ 55.0	
Duty cycle		100%		100%		100%		100%	
PDP drying performance at 7 bar		-		-		≤ +3°C		≤ +3°C	
Time to fill air tank from 0 to 9 bar	s	96		96		123		123	
Dimensions (net) w x d x h	mm	460x500x720		560x675x877		460x565x720		560x675x877	
Net weight (****)	kg	58 <sup>(***)</sup> 56 <sup>(***)</sup>	58 <sup>(***)</sup>	96 <sup>(***)</sup> 94 <sup>(***)</sup>	96 <sup>(***)</sup>	63 61	63	101 99	101
Classification under EN 60601-1		Class I.							



## Notes:

(\*) When ordering, state the version of the compressor

(\*\*) For other range of pressure consult with supplier

(\*\*\*) Weight of compressors with a KJF-1(KJF-2) unit \_add 3(3.5) kg

(\*\*\*\*) Information about the weight is for informative purposes only and applies to the product without any additional accessories.

### 7.1. Free air delivery (FAD) correction due to elevation

FAD correction table

Elevation [mamsl]	0 - 1500	1501 - 2500	2501 - 3500	3501 - 4500
FAD [l/min]	FAD x 1	FAD x 0.8	FAD x 0.71	FAD x 0.60

FAD („Free Air Delivery“) output related to conditions:

Elevation: 0 MASL

Temperature : 20°C

Atmospheric pressure : 101325 Pa

Relative humidity : 0%

## 8. PRODUCT DESCRIPTION

### 8.1. Variants

Compressor models are designed in the following variants:

**DK50 PLUS** - a base-mounted compressor for stand-alone installation in a given room

**DK50 PLUS/K** - a base-mounted compressor with a condensation and filtration unit

**DK50 PLUS/M** - a base-mounted compressor with an air dryer.

**DK50 PLUS S** - a compressor in a cabinet with effective soundproofing for in-office installations

**DK50 PLUS S/K** - a compressor in a cabinet with a condensation and filtration unit

**DK50 PLUS S/M** - a compressor in a cabinet with an air dryer

**DK50 2V** - a base-mounted compressor for stand-alone installation in a given room

**DK50 2V/K** - a base-mounted compressor with a condensation and filtration unit

**DK50 2V/M** - a base-mounted compressor with an air dryer.

**DK50 2VS** - a compressor in a cabinet with effective soundproofing for in-office installations

**DK50 2VS/K** - a compressor in a cabinet with a condensation and filtration unit

**DK50 2VS/M** - a compressor in a cabinet with an air dryer

**DK50 PLUS**

**DK50 2V**

**DK50 PLUS S  
DK50 2VS**

## 8.2. Accessories

Accessories not included in the standard order must be ordered separately!

### 8.2.1. Automatic condensate drain

The automatic condensate drain (AOK) automatically drains condensate from the compressor's air tank based on a pre-set time interval. The condensate drain (AOK) is a suitable accessory for compressor models without dryers.

Type	Use	Kit article no.
AOK 11	DK50 PLUS, DK50 2V	447000001-047

### 8.2.2. Regulator set

Compressors may be equipped with a regulator set on the compressed air outlet if specified. The regulator sets are suitable accessories for all the compressors specified above.

Type	Use	Kit article no.
REG11	DK50 PLUS, DK50 PLUS/M DK50 2V, DK50 2V/M	447000001-077

### 8.2.3. Filters set

Compressors may be equipped with a filters set on the compressed air outlet if specified. The set of filters may be equipped with a pressure regulator. Filters set are suitable accessories for all the compressors specified above.

NOTE : If a higher level of air filtration is required, this specification must be agreed on with the supplier and made clear in the order.

Type	Use	Level of filtration / $\mu\text{m}$ /	Pressure regulator	Kit article no.
FS 30FR	DK50 PLUS DK50 2V	5	yes	447000001-079
FS 30M		5 + 0.3	no	447000001-080
FS 30MR		5 + 0.3	yes	447000001-081
FS 30S		5 + 0.3 + 0.01	no	447000001-082
FS 30SR		5 + 0.3 + 0.01	yes	447000001-083
FS 31S	DK50 PLUS/M DK50 2V/M	0.3 + 0.01	no	447000001-086
FS 31SR		0.3 + 0.01	yes	447000001-087
FS 31SR2		0.01	yes	447000001-130

### 8.2.4. Condensation and filtration unit (KJF)

The compressor may be additionally fitted with a condensation and filtration unit (KJF-1, KJF-2 or KJFR-1). The KJF-1, KJF-2 or KJFR-1 ensure that the compressed air from the air tank is cooled in the cooler and the condensate is captured in the filter and automatically separated from the pneumatic distribution system. The compressed air is filtered at the same time.

Type	Use	Level of filtration / $\mu\text{m}$ /	Pressure regulator	Kit article no.
KJF-1	DK50 PLUS DK50 2V	5	no	450001011-001
KJFR-1			yes	450001011-002
KJF-2			no	450001021-001

### 8.2.5. Equipotential bonding socket

The socket allows an equipotential bonding. (Fig. 12)

Type	Use	Article no.
POAG-KBT6-EC	DK50 PLUS S, DK50 PLUS S/M	033400075-000
FLEXI-S/POAG-HK6	DK50 2VS, DK50 2VS/M	034110083-000

## 9. PRODUCT FUNCTIONALITY

### Compressor (Fig. 1)

The air pump (1) draws in atmospheric air through an inlet filter (8) and compresses it through a non-return valve (3) into an air tank (2) from which the device draws compressed air. If the pressure in the air tank drops to the switch-on pressure, the pressure switch (4) turns on the compressor and the compressor supplies compressed air into the air tank until it reaches the switch-off pressure, when the compressor switches off. The pressure hose is vented through the relief solenoid valve (13) once the air pump is shut off. The safety valve (5) prevents the pressure in the air tank from rising above the maximum allowed value. The drain valve (7) drains condensate from the air tank. Compressed, oil-free filtered air is stored in the air tank ready for use.

Condensate must be drained from the air tank at regular intervals (see Chapter 19.1).

### Compressor with membrane dryer (Fig. 3)

The air pump (1) draws in atmospheric air through an inlet filter (8) and compresses it, feeding it through the coolers (14) and the filter (15) into the dryer (9) and the dry, clean air is then fed through the non-return valve (3) into the air tank (2). A part of the air is released from the dryer with the captured moisture, which is manifested as a light air stream along the body of the dryer (9). Condensate from the filter is automatically drained at regular intervals into the collecting bottle via the condensate drain solenoid valve (16). The dryer ensures continuous drying of the compressed air. The drain valve (7) drains condensate from the air tank when drying performance is checked. Compressed, oil-free filtered air is stored in the air tank ready for use.

The pressure vessel does not need to be drained.

### Compressor with condensation and filtration unit (Fig. 2)

The air pump (1) draws in air through an inlet filter (8) and compresses it through a non-return valve (3) into the air tank (2). Compressed air from the air tank flows through a cooler (10) that cools the compressed air. The condensed moisture is trapped in the filter (11) and automatically separates as condensate (12) into the collecting bottle. Compressed, oil-free filtered air is ready for use.

Condensate must be drained from the air tank at regular intervals (see Chapter 19.1).

### Compressor cabinet (Fig. 4)

The soundproof cabinet provides compact covering of the compressor, allowing sufficient exchange of cooling air. With its design it can be placed in a dentist's office as a part of furniture. The fan under the air pump of a compressor provides cooling of compressor and it is in operation at the same time with an engine of the compressor. After prolonged use the temperature in the cabinet may rise above 40°C, causing the cooling fan blower to automatically turn on. After cooling the cabinet area down to 32°C the fan blower turns off automatically. Door of the casing with right opening may be changed to left opening (see Chapter 12.2).


**Make sure that nothing impedes the free flow of air under and around the compressor. Never cover the hot air outlet on the top back side of the cabinet.**

**If placing the compressor on a soft floor such as carpet, create space for ventilation between the base and floor or the cabinet and floor, e.g. underpin the footings with hard pads.**

**The 8-10 bar compressor model is equipped with an hour counter.**

Fig. 1 - DK50 PLUS - Compressor

Description for Figures 1-4

1. Air pump
2. Air tank
3. Non-return valve
4. Pressure switch
5. Safety valve
6. Pressure gauge
7. Drain valve
8. Inlet filter
9. Dryer
10. KJF - cooler
11. Filter
12. Condensate outlet
13. Solenoid valve
14. Dryer cooler
15. Filter
16. Condensate drain solenoid valve
17. Compressed air outlet
18. Circuit breaker switch
19. Bottle
20. Magnetic holder
21. Fan
22. Cabinet fan
23. Door pin
24. Rectifying screw
25. Lock
26. Supporting bracket
27. Wall stop
28. Switch
29. Pressure gauge
30. Door hinge
31. Handle
32. Casters
33. Pressure gauge hose
34. Socket
35. Cabinet connector with cord

Fig. 2 – DK50 PLUS/K - Compressor with condensation and filtration unit KJF-1


Fig. 3 - DK50 2V/M - Compressor with dryer




Fig. 4 - Cabinet


## 10. PNEUMATIC DIAGRAM

DK50 PLUS, DK50 PLUS S, DK50 2V, DK50 2VS

DK50 PLUS/M, DK50 PLUS S/M, DK50 2V/M, DK50 2VS/M

### Legend for compressed air diagram

- |                     |                                     |
|---------------------|-------------------------------------|
| 1. Inlet filter     | 9. Safety valve                     |
| 2. Air pump motor   | 10. Air tank                        |
| 3. Fan              | 11. Drain valve                     |
| 4. Relief valve     | 12. Cooler                          |
| 5. Noise silencer   | 13. Coalescence filter              |
| 6. Non-return valve | 14. Membrane dryer                  |
| 7. Pressure gauge   | 15. Condensate drain bottle         |
| 8. Pressure switch  | 16. Condensate drain solenoid valve |
|                     | 17. Outlet valve                    |

## **INSTALLATION**

### **11. CONDITIONS FOR USE**

- Compressors may only be installed and operating in dry, well-ventilated and clean environments with environmental parameters that meet the requirements specified in Chapter 7, Technical data. The compressor must be installed so that it is accessible at all times for operating and maintenance. Please ensure that the label on the device is readily accessible.
- The compressor must stand on a flat, sufficiently stable base (be aware of the weight of the compressor, see Chapter 7, Technical data).
- Compressors may not be operated outdoors or in otherwise wet or damp environments. Do not use the compressor in the presence of explosive gases, dust or flammable liquids.
- Before connecting the compressor to medical equipment, the supplier must confirm that it meets all requirements for its use. Refer to the technical data of the product for this purpose. When a unit is to be built-in, classification and evaluation of compatibility must be done by the manufacturer or supplier of the product to be used.
- Any other use or use outside this framework is not considered as intended use. The manufacturer is not responsible for any damages arising from such use. The operator/user assumes all risks.

**Only a qualified professional may install the compressor and place it into operation for the first time. This professional is obliged to train operating staff as to the use and maintenance of the device. Installation and training of all operators shall be confirmed by the installer's signature on the certificate of installation.**

**Prior to installation, ensure that the compressor is free of all transport packaging and stabilizers to avoid any risk of damage to the product.**

**Parts of the air pump may be hot and reach hazardous temperatures during compressor operation and may pose a contact risk for operators or materials. Burn or fire hazard! Caution! Hot surface!**

### **Ambient operating conditions**

*Temperature:* +5°C to +40°C,  
*Max. relative humidity:* 70%,  
*Max. absolute humidity:* 15 g/m<sup>3</sup>.

## 12. PLACEMENT OF THE COMPRESSOR

The device shall only be installed by a qualified professional.

- Unwrap the compressor from the packaging.

### 12.1. Handling and releasing the compressor

- Position the compressor at the site of future operation (Fig. 5)


Fig. 5

- Remove the transport securing elements/bolts from the air pumps. (Fig. 6)

**Remove all devices used to secure the air pumps once the compressor set has been installed and mounted at the site of final installation!**



Fig. 6

## 12.2. Changing door direction

- Remove the door, rectification screw (3) and bracket (2) for hinge D (4).
- Install the bracket for hinge D (2) on the left side of the cabinet.
- Rotate the door 180°.
- Insert the spacer (5) between hinge H (1) and the bottom of the door.
- Install the door.
- Remove the door lock (6) and rotate 180°.
- Remove the latch (7) and rotate 180°.
- Install the lock.


### 12.3. Adding foam to the cabinet


Cabinet-mounted compressors **not equipped with a dryer** require foam be added to the tunnel in the cabinet.

This foam is provided in the basic equipment delivered with the product.


A

B

C

## 13. PNEUMATIC CONNECTIONS

### 13.1. Compressed air outlet (Fig. 7)

Connect the pressure hose (2) to the compressed air outlet fitting (1) on the compressor. Run the hose to the compressed air manifold or directly to the dental set.


Fig. 7

Route the pressure hose through the opening in the rear wall of the cabinet for cabinet-mounted compressors. (Fig. 9)

### 13.2. Connecting the cabinet pressure gauge hose to the compressor (Fig. 8)

Remove the threaded plug from the pneumatic block on the compressor. Connect the cabinet pressure gauge hose to the threaded fitting.


Fig. 8

### 13.3. Condensate outlet (Fig. 9)

- Connect the condensate drain hose to the condensate collecting bottle on compressors with dryers.

Fig. 9

- Route the hose through the opening in the rear wall of the cabinet for compressor models with a dryer.

### 13.4. Floor installation

(Fig. 10)

Connect the compressor using the hoses prepared in advance in the floor, as per the installation plan.


Fig. 10

### 14. ELECTRICAL CONNECTIONS

The product is equipped with a grounded plug. Insert the mains plug into a rated mains socket.

- Keep the socket easily accessible to ensure that in an emergency the device can be safely disconnected from the mains.
- Connect the cabinet to the compressor using the electrical cord with the connector (1) into the socket (2). (Fig. 11)
- Route the power cord through the opening in the rear wall of the cabinet for cabinet-mounted compressors.


Fig. 11

**Ensure full compliance with all local electrical codes. The mains voltage and frequency must comply with the data stated on the device's label.**

- Connection to the power distribution cabinet must be max.16 A.
- Connect the pin for equipotential bonding  $\varnothing$  6 mm (1) to the electrical circuit using the defined method under valid electrical engineering regulations. The socket for equipotential bonding (2) is an accessory and is not provided in the basic package. (Fig. 12)


Fig. 12

**Ensure the electrical cable does not touch hot compressor components. Risk of electric shock!**

**If any electrical cord or air hose is damaged it must be replaced immediately.**


### 14.1. Jumper removal

Remove the jumper from the terminal strip in the electrical panel for cabinet-mounted compressors. **The switch on the compressor cabinet will not work properly if the jumper is not removed!**


- Remove the cover on the electrical panel.
- Remove the jumper from the terminal strip.
- Reinstall the cover on the electrical panel.

Keep the removed jumper for service purposes (see Chapter 19.11 - Jumper connections).

Compressor with jumper installed


230V


400V

Compressor with jumper removed

230V

400V


## 15. CONNECTION DIAGRAM

**DK50 PLUS, DK50 PLUS S, DK50 2V, DK50 2VS** 5-7 bar, 6-8 bar, 8-10 bar

1/N/PE ~ 230V 50 Hz  
~ 230V 60 Hz  
~ 115V 60 Hz  
ELECTRIC OBJECT OF 1st. CAT.

NOTE:

\* - only 8-10 bar


**DK50 PLUS/M, DK50 PLUS S/M, DK50 2V/M, DK50 2VS/M**      **5-7bar, 6-8 bar, 8-10 bar**

1/N/PE ~ 230V 50 Hz  
~ 230V 60 Hz  
~ 115V 60 Hz  
ELECTRIC OBJECT OF 1st. CAT.

NOTE:

\* - only 8-10 bar

**DK50 2V, DK50 2VS      5-7 bar, 6-8 bar, 8-10 bar**


1/N/PE ~ 400V 50Hz

MAINS TN-S [TN-C-S]

ELECTRIC OBJECT OF 1st. CAT.

NOTE:

\* - only 8-10 bar

**DK50 2V/M, DK50 2VS/M     5-7bar, 6-8 bar, 8-10 bar**


1/N/PE ~ 400V 50Hz

MAINS TN-S [TN-C-S]

ELECTRIC OBJECT OF 1st. CAT.


NOTE:

\* - only 8-10 bar

## Compressor cabinet

1/N/PE ~ 230V 50Hz  
~ 115V 60Hz  
ELECTRIC OBJECT OF 1st. CAT.

### Legend for electrical diagram

M1	Compressor motor	C1	Capacitor
E1	Compressor fan	B1	Pressure switch
E2	Dryer fan	X1	Terminal box
M10	Relief valve	F1	Breaker switch
B2	Temperature switch	M4	Condensate drain valve
E10	Cabinet fan	P21*	Hour counter
X10, X2	Connector	S10	Switch
F11	Breaker	P11	Indicator
Q11	Contact	X2	Socket

Note:

\* - The 8-10 bar compressor model is equipped with an hour counter.

**\*\* - Only connect the jumper for compressors without cabinet (Chapter.19.11)**

**OPERATION**

**IN CASE OF EMERGENCY, DISCONNECT THE COMPRESSOR FROM THE MAINS (PULL OUT THE MAINS PLUG).**

**THE AIR PUMP HAS HOT SURFACES.  
BURNS OR FIRE MAY RESULT IF CONTACT IS MADE.**

**During prolonged operation of the compressor, the temperature in the cabinet may increase to over 40°C. At this point the cooling fan automatically switches on. After cooling the space to under 32°C, the fan switches off.**

**Automatic start: when pressure in the pressure tank decreases below the switch-on pressure, the compressor automatically switches on. The compressor automatically switches off when pressure in the air tank reaches the switch-off pressure.**

**The working pressure settings for the pressure switch set by the manufacturer cannot be changed. Compressor operation at a working pressure below the switch-on pressure indicates high air consumption (see the chapter Malfunctions)**

**Required drying performance can only be achieved when following the defined operating conditions!**

**Drying performance will decline and the achieved dew point will drop if the dryer is operated at any pressure below the minimum working pressure!**

**WHEN THE DRYER IS OPERATED AT AMBIENT TEMPERATURE THAT IS HIGHER THAN MAXIMUM OPERATING TEMPERATURE, THE DRYER MAY BE DAMAGED!**

**16. COMMISSIONING**

- Make sure all transport stabilizers were removed.
- Check that all pressurized air line connections are secure.
- Check to ensure power is connected correctly to the compressor.
- Check the position of the breaker switch to ensure it is in the "I" position. Turn the breaker switch (4) to the "I" position if found in the "0" position. (Fig. 13)
- For cabinet-mounted compressors, turn the switch (5) on the front side of the cabinet to the "I" position; a green indicator indicates the device is in operation. (Fig. 13)
- Check the connection of the cabinet connector to the compressor. (Fig. 11)
- Check to ensure the jumper has been removed from the terminal strip in the electrical panel for cabinet-mounted compressors. (chap.14.1.)
- Check the connection of the the cabinet pressure gauge hose to the compressor pneumatic block. (Fig. 4 and 8)

**The compressor is not equipped with a backup power supply.**

## 17. SWITCHING THE COMPRESSOR ON

(Fig. 13)

Start the compressor at the pressure switch (1) by turning the switch (2) to the "I" position. Check the position of the breaker switch to ensure it is in the "I" position. If not, turn the breaker switch (4) to the "I" position. Turn on the switch (5) on the front of the cabinet for compressors in the cabinet, the indicator will turn green. The compressor begins running, the air tank is filled, and the switching pressure and pressure switch turn the compressor off. The compressor operates in automatic mode, switched on and off by the pressure switch, depending on compressed air consumption. Check the values of the switch-on and switch-off pressure using the pressure gauge (3). A tolerance of  $\pm 10\%$  is acceptable. The pressure in the air tank must not exceed the permitted operating pressure.


Fig. 13

**Adjustments to the pressure range on the compressor's pressure switch are prohibited. The pressure switch (1) was set up by the manufacturer and only a qualified technician trained by the manufacturer may make any changes to its settings.**

**Compressor** – when first started and placed into service, the compressor fills the air tank until the switching pressure and the pressure switch turn off the compressor. The compressor operates in automatic mode, switched on and off by the pressure switch, depending on compressed air usage.

**Compressor with dryer** – the compressor operates in the same manner as above, only the compressed air passes through a cooler and a dryer that removes moisture from the compressed air.

**Compressor with a condensation and filtration unit (KJF)** - during usage by the device, the compressed air passes through the KJF unit, where the air is cooled, filtered and condensed liquid is captured and automatically drained into a collecting bottle.

## 18. COMPRESSOR SHUT-DOWN

(Fig. 13)

Compressor shut-down for either service or any other reason is performed at the pressure switch (1) by turning the alternation switch (2) to the "0" position and **pulling the mains plug from the socket**. This disconnects the compressor from power supply. Then open the drain valve (Fig. 14) to decrease the pressure in the air tank to zero.

## MAINTENANCE

### 19. DEVICE MAINTENANCE

#### **Warning!**

The operator shall ensure completion of repeated testing of the device at least once every 24 months (EN 62353) or at intervals defined by applicable national legal regulations. A record of these test results shall be completed (e.g. per EN 62353, Annex G) together with the measurement methods.

The device has been designed and manufactured to keep maintenance to a minimum. The following work must be performed to retain the proper and reliable operation of the compressor.

**Before starting compressor maintenance work, it is necessary to check if the compressor can be disconnected from the device to ensure that the person using the device is not at risk of health damage and there is no risk of any other material damages!**

**Air pump components (head, cylinder, pressure hose, etc.) are very hot during and shortly after compressor operation – do not touch these components!**

**Repair works beyond normal maintenance can be performed only by qualified personnel or the manufacturer's customer service.  
Use only spareparts and accessories approved by the manufacturer.**

**PROTECT EYESIGHT, WEAR GOGGLES, WHEN VENTING COMPRESSED AIR FROM THE COMPRESSED AIR CIRCUIT (AIR TANK).**

**The 8-10 bar compressor model is equipped with an hour counter.**

The work below may only be performed by trained personnel as follows:

**BEFORE STARTING ANY MAINTENANCE WORK, FIRST DISCONNECT THE COMPRESSOR FROM THE MAINS (PULL OUT THE MAINS PLUG) AND VENT THE COMPRESSED AIR FROM THE AIR TANK.**

**Allow the equipment to cool down before maintenance, service or connecting/disconnecting the compressed air supply!**

**Reconnect the grounding lead if removed during service work to its previous location once the work is complete.**

## 19.1. Maintenance intervals

Time interval ***	once a day	once a week	once a year	once every 2 years	2000 hours	4000 hours	6000 hours	8000 hours	10000 hours	12000 hours	Chapter	Set of spare parts	Performed by
Product operational checks	x										19.2	-	user
Drain the condensate from the air tank - At high RH )** - Pour the condensate out of the vessel )****	x										19.5	-	user
Drain the condensate from the air tank - At normal RH )**		x									19.5	-	user
Product function check		x									9	-	user
Leak check connections and inspect device			x								19.3	-	qualified professional
Inspection of electrical connections			x								19.4	-	qualified professional
Cooler and fan check			x								19.10	-	qualified professional
Filter element replacement in the dryer filter			x								19.8	025200304-000	qualified professional
Filter element replacement in the KJF			x								19.9	025200061-000	qualified professional
Safety valve check			x								19.6	-	qualified professional
Conduct a "repeated test" per EN 62353				x							19	-	qualified professional
Replacing the inlet filter on the DK50 PLUS air pump )*				x		x		x		x	19.7	025200126-000	qualified professional
Replacing the inlet filter and pre-filter on the DK50 2V air pump )*				x		x		x		x	19.7	025200139-000 025200150-000	qualified professional

)\* data is presented in hours; if not possible, then the data is considered in years

)\*\* only for compressors without dryers

)\*\*\* For 60 Hz compressor models, reduce the time interval in hours by 20%

(2000 h/1600 h, 4000 h/3200 h, 6000 h/4800 h, 8000 h/6400 h, 10000 h/8000 h, 12000 h/9600 h)

)\*\*\*\* only for compressors with dryers

### **19.2. Operational checks**

- Check air pump condition – the air pump s should be operating normally without excessive vibration or noise. Troubleshoot any problem or call in service personnel if trouble is detected
- Visually inspect fan operation – the fans must be operating when the air pumps are running. Troubleshoot any problem or call in service personnel if trouble is detected
- Check to ensure the power cable and pneumatic hoses are undamaged. Replace damaged components or call in service personnel.
- Check the ambient temperature – the ambient temperature must be below the temperature limit (40 °C). Cool the space if the temperature is high.
- On the compressor with dryer – open the plug on the condensate collecting bottle to drain the condensate.

### **19.3. Pneumatic connection leakage check and device inspection**

#### **Leakage check:**

- Perform the leakage check of pneumatic connections when the compressor is operating – pressurizing.
- Use a leakage analyser or soapy water to check all joints and connections for leaks. Tighten or reseal the connection where leaks are found.

#### **Device inspection:**

- Check the condition of the air pump for normal operation and noise levels.
  - Fan operation check - the fans must be running during the defined compressor work cycles
  - Check of temperature switch operation (B2) – heat the temperature switch to a temperature above 40°C (e.g. using a heat gun, being careful not to warp any plastic pieces nearby). Fan EV1 (and EV2 for a compressor with dryer) will start up once the temperature reaches 40°C so long as power is connected to the compressor.
  - Check filter condition – the filters must be undamaged and sufficiently clean.
  - Check the condition of the air pump itself - ensure there is no contamination inside the crankcase or play in the crankshaft
- Replace any defective parts as needed.

### **19.4. Inspection of electrical connections**

**Check all electrical connections on the device with the mains disconnected!**

#### **Inspection**

- Check the mechanical function of the main switch.
- Check to ensure the power cable and connected wires are undamaged.
- Visually inspect the connection of individual cables to the terminal strip.
- Inspect all screw terminals for the protective green and yellow PE grounding conductor.

### 19.5. Condensate drain

#### Compressors (Fig. 14)

Draining the condensate from the pressure vessel is recommended during regular operation. Switch off the compressor at the mains. Reduce air pressure in the device to max. 1 bar by releasing air via a connected device. Install the hose from the basic package into the drain valve. Place a vessel under the hose and open the drain valve (1) and drain the condensate from the tank.


Fig. 14

Carefully and slowly open the drain valve. Quickly opening the drain valve is accompanied by a loud sound and uncontrollable stream of accumulated condensate.

#### Compressors with condensation and filtration unit

(Fig. 15)

During regular use, condensation is automatically drained via the drain valve of the condensation unit filter. To check that the automatic drain is working properly, open the valve (2) of the drain vessel (1) by turning to the left. Release a small amount of condensate from the vessel. Close the valve (2) by turning to the right.


Fig. 15

Both compressor types can be equipped with an automatic condensate drain on the air tank to allow condensate to drain without operator intervention (see the— Accessories chapter).

#### Compressors with air dryer (Fig. 16)

Condensate from compressors with air dryers is automatically drained into a collecting bottle. The bottle must be drained regularly.

Fig. 16

Monitor the level in the vessel using the 1 l or 2 l markings (depending on the volume of the vessel) at least once a day. The vessel may overflow if the condensate is not drained out of the vessel at the defined interval.

**THE FOLLOWING IS REQUIRED BEFORE ANY FOLLOW-UP CHECKS:**

**For cabinet-mounted compressors,** unlock the door and open the door on the cabinet.

### 19.6. Safety valve check

(Fig. 17)

When the compressor is operated for the first time, make sure that the safety valve is working properly. Turn the screw (2) on the safety valve (1) several rotations to the left until the safety valve releases air. Let the safety valve vent for only a few seconds. Turn the screw (2) to the right until it seats, closing the valve.


Fig. 17

Never use the safety valve to vent air pressure in the air tank. It could damage the safety valve. The valve is set to the maximum permitted pressure by the manufacturer. Adjustments are not permitted!

**WARNING! COMPRESSED AIR CAN BE DANGEROUS. WEAR EYE PROTECTION WHEN VENTING ANY AIR. COMPRESSED AIR MAY DAMAGE THE EYES OR CAUSE VISION PROBLEMS.**

### 19.7. Replacing the inlet filter and pre-filter

(Fig. 18)

The inlet filter is located inside the compressor crankcase cover.

Inlet filter replacement:

- Pull out the rubber plug by hand (2).
- Remove the used filter (1).
- Insert a new filter and replace the rubber plug.

Pre-filter replacement:

- Pull out the pre-filter by hand (3).
- Replace the element and re-insert.



Fig. 18

### 19.8. Filter element replacement

(Fig. 19)

- Remove the hose (1) from the quick connector.
- Use a wrench (2) to release the filter vessel (3) and remove.
- Pull down on the filter element (4) to remove.
- Insert a new filter element.
- Re-install the filter vessel.
- Gently tighten the filter vessel with the wrench.
- Re-install the hose on the quick connector.




Fig. 19

### 19.9. Filter element replacement in the KJF unit

**Before proceeding, vent the air tank to zero pressure and disconnect the device from power.**

(Fig. 20)

- Release the catch (1) on the filter by pulling down, rotate the vessel (2) and pull out.
- Pull out the holder with filter (3), rotate and take out from the vessel.
- Rotate the filter retainer (4).
- Replace the filter element (5) and reinstall the filter retainer (4), secure by rotating.
- Place the filter bracket (3) back into the vessel and rotate to secure it in the place.
- Reinstall the filter housing and rotate to secure it until the catch clicks.

Fig. 20

### 19.10. Cooler and fan check (Fig.3)


The equipment altogether, and especially the compressor fan, cooler fan (21) and cooler (14), must be kept clean to ensure efficient drying performance; vacuum or blow down the cooling fins and fans with compressed air to remove any dust from the surface.

### 19.11. Jumper connections


It may be necessary to start the cabinet-mounted compressors outside of the cabinet for service purposes. This means starting the compressor without the cabinet connector attached. To start the compressor, a jumper must be installed in the electrical box.

- DK50 PLUS S and DK50 2VS compressor models require **the electrical jumper be installed** to start the compressor outside the cabinet. (Fig. 21)

Fig. 21

230V


400V


- DK50 PLUS S and DK50 2VS compressor models require **the electrical jumper be removed** to start the compressor inside the cabinet. (Fig. 22)

Fig. 22

230V

400V

## TROUBLESHOOTING

**CAUTION! BEFORE PROCEEDING, DEPRESSURIZE THE AIR TANK TO ZERO AND DISCONNECT THE DEVICE FROM THE MAINS.**

For permanently high efficiency of drying, it is necessary to maintain the whole device, and mainly fan clean – regularly clean the surface of fan and cooling fins of cooler.

Troubleshooting can be performed only by qualified personnel.

MALFUNCTIONS	POSSIBLE CAUSE	REMEDY
Compressor does not start	No voltage at the pressure switch  Motor winding failure, damaged thermal overload protection Faulty capacitor Seized piston or other rotating part Pressure switch does not switch	Check voltage in socket Check the breaker switch - and switch to the "I" on position Loose conductor at terminal – tighten Check the electrical cord – replace defective cord Replace motor or windings  Replace capacitor Replace damaged parts Check the function of the pressure switch
Compressor often switches on	Air leak in compressed air distribution system Leaking non-return valve  Large volume of condensed liquid in pressure vessel Low compressor performance	Check compressed air distribution system – seal loose joints Clean non-return valve, replace seals, replace non-return valve Drain condensed liquid Check the time to fill the air tank

Low pressure in air tank (compressor running constantly)	High air usage by device, Leak from pneumatic distribution system, Low air pump output Air pump malfunction Dryer malfunction	
Prolonged operation of the compressor	Leak from pneumatic distribution system Worn piston ring Inlet filter is plugged Defective solenoid valve	Check compressed air distribution system – seal loose joint Replace worn piston ring Replace old filter with a new filter Repair or change the valve or coil
Compressor is noisy (knocking, metal noises)	Damaged piston bearing, piston rod, motor bearing Loose or cracked spring	Replace damaged bearing Replace damaged spring
Dryer doesn't dry (condensed water in the air) *	Cooling fan not working	Replace the fan Check the power source
	Damaged dryer	Replace the dryer
	Non-functioning automatic condensate drain	Clean / replace

)\* After a dryer malfunction, the interior surfaces of the air tank must be thoroughly cleaned and all condensed liquid removed.

**Check the dew point of the air leaving the air tank (see Chapter 7 - Technical data) to protect the device from damage!**

## 20. REPAIR SERVICE

Guaranteed and post-guarantee repairs must be done by the manufacturer, its authorized representative, or service personnel approved by the supplier.

**The manufacturer reserves the right to make changes to the device without notice. Any changes made will not affect the functional properties of the device.**

## 21. STORAGE

If the compressor will not be used for a prolonged time period, drain all condensate from the air tank. Then turn on the compressor for 10 minutes, keeping the drain valve open (1) (Fig. 14). Switch off the compressor using the switch (2) at the pressure switch (1) (Fig. 13), close the drain valve and disconnect the device from the mains.

## 22. DISPOSAL OF THE DEVICE

- Disconnect the device from the mains.
- Release air pressure in the pressure tank by opening the drain valve (1) (Fig. 14).
- Dispose of the device following all applicable environmental regulations.
- Entrust a specialised company to sort and dispose of waste.
- Worn out components have no negative environmental impact.

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